

Introducing DermaFlush™

An economic, convenient and easy-to-travel form of saline when added to water.

INGREDIENTS:

Premixed Saline Packets: Sodium Chloride, Boric Acid.

- Soothing
- Preservative free
- Iodine free
- Isotonic solution or hypertonic solution
- Reduces edema
- Draws out mucous and debris
- For small, exotic and large animal practice



INDICATIONS: Use anywhere a saline flush may be beneficial, including ears, sinuses, wounds and external surgery.

Wound Flush: If using the Earigator, clean and flush the unit and the tubing with a mild bleach solution and water prior to use. Cut the tip off the catheter tipped end of a 10 - 12 Fr red rubber catheter to fit the trumpet valve nipple. Then fill the Earigator flush bottle with a saline solution made with DermaFlush. The catheter can be used to snake under the skin and into loose skin or muscle tissue or abscesses and fistulated areas with fluid under pressure provided by the Earigator.

Ear Flush: DermaFlush can be used to flush ears. It can be used in the Earigator to flush ears while the patient is under sedation. It is gentle and non-irritating.

Nasal Flush: Nasal flushing may be helpful in cases of chronic nasal discharge. Sedate the patient and flush the sinus area with hypertonic saline. This helps reduce edema and gently removes debris and discharge from the sinuses.

Available in boxes of single use packets.



www.dermazoo.com / info@dermazoo.com

Using the Earigator

To start, make sure the unit is powered OFF. Unscrew the irrigation bottle and fill it with irrigation solution to 800 mL maximum. The bottle should then be replaced in the pocket, screwing it firmly into the cap, and the desired irrigation tip should be installed on the hand piece fitting. The unit may then be powered ON.



Adjust the irrigation pressure to between 5 and 15 PSI by pulling out the regulator knob or pulling out the red locking ring (depending on your model), and turning it clockwise for more pressure or counter-clockwise for less. Note that the regulator should be reduced below the set pressure and increased slowly to achieve the desired setting.

Press and hold the BLUE button to discharge irrigation solution. Press and hold the WHITE button to suction solution back. You may adjust irrigation and suction pressures to suit your tip size and preferences.